



PARKER MEDICAL SUPPLIES



DESTANDAU ENDOSPINE SET

Product Catalog



+1 (832) 343-3140



info@parkermedicalsupsies.com



www.parkermedicalsupsies.com

2800 Post Oak Blvd, Ste 5600, Houston TX 77056, United States






PMS0040

TELESCOPE

Straight Forward Telescope 0°, enlarged view, diameter 4 mm, length 18 cm, autoclavable, fiberoptic light transmission incorporated



PMS0041

LOCALIZATION DEVICE

for fluoroscopic determination of the point of incision for ENDOSPINE® Operating Tube 28163 DW



PMS0042

BONE PUNCH

dismantling, 90° upbiting, not through-cutting, 3 mm, working length 18 cm



PMS0043

BONE PUNCH

dismantling, upbiting, 45° forward, not through-cutting, 3 mm, working length 18 cm



PMS0044

TREPINE

with round handle, diameter 3 mm, working length 22 cm



SPOON FORCEPS

dismantling, robust, oval, single action jaws, spoon size 3 x 10 mm, working length 15 cm

BIPOLAR COAGULATION FORCEPS

with connecting pin for bipolar coagulation, width of jaws 1 mm, diameter 5 mm, working length 20 cm,



PMS0045

CHISEL

flat, straight, with handle, distal width 15 mm, working length 9 cm



PMS0046

ELEVATOR

spatula slightly curved, distal width 5 mm, working length 13 cm



PMS0047

PALPATION HOOK

blunt, distally angled 90°, hook length 5.5 mm, working length 13 cm



PMS0048

SUCTION TUBE

angled, with cut-off hole, diameter 3.7 mm, working length 11 cm



PMS0049

BIPOLAR HIGH FREQUENCY CORD



PMS0050



DESTANDAU OPERATING TUBE WITH WORKING INSERT

Operating Tube, oval, with Obturator for use with Working Insert with positioning detent, with adjustable Nerve Protector for use with Operating Tube with working channel diameter 8 mm and irrigation channel

PMS0051



PARKER MEDICAL SUPPLIES

CONTACT US



+1 (832) 343-3140



info@parkermedicalsupplies.com



www.parkermedicalsupplies.com

VISIT US:

2800 Post Oak Blvd, Ste 5600, Houston TX 77056, United States