



PARKER MEDICAL SUPPLIES



PMS HYSTROSCOPY AND UROLOGY INSTRUMENTS

Product Catalog



+1 (832) 343-3140



info@parkermedicalsupplies.com




www.parkermedicalsupplies.com

2800 Post Oak Blvd, Ste 5600, Houston TX 77056, United States

HYSTEROSCOPY & UROLOGY

WORKING ELEMENT MONOPOLAR & BIPOLAR (COMBO)

- Monopolar Resectoscope can be used in complex intrauterine operations without any difficulty
- Bipolar Resectoscope are particularly well suited for infertility patients
- Bipolar resection benefits from generation of plasma to achieve outstanding cutting power
- By using HF current several operations like hysteroscopic endometrium ablation, Myoma resection, adhesiolysis or septum dissection can be performed
- An electrolyte free irrigation fluid is used as distension media



Monopolar / Bipolar Passive Resectoscope Complete (For TCRE/TURP) **PMSMBPRC-4/302**



Outer sheath




Inner Sheath



Electrode Cutting Loop Monopolar/ Bipolar



Working Element Active & Passive Resectoscope Monopolar / Bipolar



Resectoscope Sheath, Continuous flow with rotatable inner sheath and blunt Obturator



Obturator for resectoscope sheath

ELECTRODES

HF-Electrodes, monopolar, single stem (for 24 fr)



Electrode Cutting Loop,



Electrode Roller Ball




Electrode KNIFE, Angled



Electrode Cutting Loop Bipolar

DIAGNOSTIC AND OPERATIVE SHEATH- CIRCULAR



Hysteroscopic Examination Sheath with Atraumatic Distal End Rotatable Stop Cock




Hysteroscopic Operative Sheath Round Continuous Flow, one working channel upto 6fr & detachable Suction irrigation plugs



Hysteroscopic Diagnostics Sheath, contributions with 2 detachable Suction Irrigation plugs for

OPERATIVE SHEATH-OVAL (DETACHABLE AND NON-DETACHABLE)



Continuous flow with Atraumatic distal end two stop cock and one working channel of (5fr) semi flexible and flexible instruments (5fr)



SEMI FLEXIBLE INSTRUMENTS



CORNUAL CANNULATION SET PMSCS-5FR



UROLOGY

STONE PUNCH





PARKER MEDICAL SUPPLIES

CONTACT US



+1 (832) 343-3140



info@parkermedicalsupplies.com



www.parkermedicalsupplies.com

VISIT US:

2800 Post Oak Blvd, Ste 5600, Houston TX 77056, United States