



PARKER MEDICAL SUPPLIES



PMS DISPOSABLE INSTRUMENTS

Product Catalog



PARKER MEDICAL SUPPLIES



Disposable Laparoscopic Instruments

Laparoscopic handheld instruments consist of linkless, stainless steel drive

mechanism which delivers more precise "hand-to-hand" operation



Features and Benefits



Large jaw aperture for greater flexibility



Ergonomic design for superior feel


























Easily maneuverable, 360° degree rotation knob



Extended insulation to prevent stray currents

Disposable Laparoscopic Instruments

Scissors (Purple Knob)	Tip Style	REF.NO.	Description
		PMSLES 1450	Curved Metzenbaum, Extended Insulation, 5mm x 330mm
		PMSLES 1460	Curved Metzenbaum, Extended Insulation, 5mm x 430mm
Graspers (Green knob, with Ratchet)		PMSLEG 1510	Maryland, 5mm x 330mm
		PMSLEG 1520	Fenestrated (Duckbill), 5mm x 330mm
		PMSLEG 1530	Fenestrated (Long), 5mm x 330mm
		PMSLEG 1540	Fenestrated, 5mm x 330mm
		PMSLEG 1550	Tapered (Dolphin), 5mm x 330mm
		PMSLEG 1560	Babcock, 5mm x 330mm
		PMSLEG 1570	Fenestrated Babcock, 5mm x 330mm
		PMSLEG 1580	Meeker Grasper, 5mm x 330mm
		PMSLEG 1590	Allis Grasper, 5mm x 330mm
	Dissectors (Blue knob, Non-ratcheting)		PMSLED 1510
		PMSLED 1520	Fenestrated (Duckbill), 5mm x 330mm
		PMSLED 1530	Fenestrated (Long), 5mm x 330mm
		PMSLED 1540	Fenestrated, 5mm x 330mm
		PMSLED 1550	Tapered (Dolphin), 5mm x 330mm
		PMSLED 1560	Babcock Grasper, 5mm x 330mm
		PMSLED 1570	Fenestrated Babcock, 5mm x 330mm
		PMSLED 1580	Meeker Grasper, 5mm x 330mm
		PMSLED 1590	Allis Grasper, 5mm x 330mm



Disposable Retrieval Bags

Multi-size, Reinforced, Leak Proof

Disposable Retrieval Bags

Disposable Retrieval Bags is one of the most economical retrieval system available in the current laparoscopy market.

The product with function of deployed automatically, easy to remove and unload during procedures.

Features and Benefits



Auto retrieval with different sizes




TPU bags for superior durability



No breaks and leaks



Exceptional safety and security



Product	Description	Packaging
PMSEB-0100	100ml, 10mm x 330mm, single use, sterile	1/pk,10/bx,100/ctn
PMSEB-0200	200ml, 10mm x 330mm, single use, sterile	1/pk,10/bx,100/ctn
PMSEB-0400	400ml, 10mm x 330mm, single use, sterile	1/pk,10/bx,100/ctn
PMSEB-0700	700ml, 10mm x 330mm, single use, sterile	1/pk,10/bx,100/ctn
PMSEB-1200	1200ml, 10mm x 330mm, single use, sterile	1/pk,10/bx,100/ctn
PMSEB-1600	1600ml, 12mm x 330mm, single use, sterile	1/pk,10/bx,100/ctn



Disposable Specimen Pouches

Simple, Cost-effective and Reliable

Disposable Specimen Pouches

The Disposable Specimen Pouch is a simple and low-cost specimen retrieval system with superior durability.

Our pouches offer easy and safe specimen capture and removal during surgical procedures.

Features and Benefits



A variety of sizes to meet different demands



Drawstring for intra-abdominal closure



Special knot to prevent re-opening



Transparent TPU bags



Exceptional safety and security



Product	Description	Packaging
PMSEP-0100	100ml, 110mm x 150mm, single use, sterile	1/pk, 100/bx, 1000/ctn
PMSEP-0200	200ml, 115mm x 152mm, single use, sterile	1/pk, 100/bx, 1000/ctn
PMSEP-0400	400ml, 160mm x 140mm, single use, sterile	1/pk, 100/bx, 1000/ctn
PMSEP-0700	700ml, 170mm x 200mm, single use, sterile	1/pk, 100/bx, 1000/ctn
PMSEP-1600	1600ml, 170mm x 290mm, single use, sterile	1/pk, 100/bx, 1000/ctn
PMSEP-3000	3000ml, 210mm x 300mm, single use, sterile	1/pk, 50/bx, 500/ctn



Disposable Retrieval Device with Memory Wire

Flexible Body, Multi-size, Anchor-hold

Disposable Retrieval Device with Memory Wire

The Disposable Retrieval Device with Memory Wire is an unique, self-opening specimen retrieval system with superior durability.

Our retrieval device offer easy and safe capture and removal during surgical procedures.

Features and Benefits



Flexible memory wire specimen bag and introducer



A variety of sizes to meet different demands



Transparent TPU bags



Exceptional safety and security



Product	Description	Packaging
PMSSEB-0100	100ml, 110mm x 150mm, 10mm Introducer, single use, sterile	1/pk,10/bx,100/ctn
PMSSEB-0200	200ml, 100mm x 130mm, 10mm Introducer, single use, sterile	1/pk,10/bx,100/ctn
PMSSEB-0400	400ml, 160mm x 140mm, 10mm Introducer, single use, sterile	1/pk,10/bx,100/ctn
PMSSEB-0700	700ml, 170mm x 200mm, 12mm Introducer, single use, sterile	1/pk,10/bx,100/ctn



Disposable Redeployable Ripstop Retrieval Bags

Extra Reinforced, Excellent Encapsulated, Leaf Proof

Disposable Redeployable Ripstop Retrieval Bags

Disposable Redeployable Ripstop Retrieval Bag is made of nylon with a thermoplastic polyurethane (TPU) coating, with characteristic of tear-resistant, impervious to fluids and multiple specimens retrieval. The bags offer efficient and safe tissue removal in surgical procedures.

Features and Benefits



Redeployable bag allows for retrieval of multiple specimens in one sugry



Closure structure prevents bag from re-opening



Pre-loaded spring automatically opens bag upon deployment



Radiopaque thread allows bag to be visible in x-rays



ReRipstop Nylon with Polymer coating for leakproof performance



Product	Description	Packaging
PMSNEB-0100	100ml,107mm x 146mm,Single Use,Sterile	1/pk, 10/bx, 100/ctn
PMSNEB-0200	200ml,120mm x 170mm,Single Use,Sterile	1/pk, 10/bx, 100/ctn



Disposable Ripstop Retrieval Bags

Extra Reinforced, Excellent Encapsulated, Leak Proof

Disposable Ripstop Retrieval Bags

Disposable Ripstop Retrieval Bag is made of nylon with a thermoplastic polyurethane (TPU) coating, with characteristic of tear-resistant and impervious to fluids.

The bags offer easy and safe tissue removal in surgical procedures.

Features and Benefits



Nylon with polymer coating for leak-proof performance



A variety of sizes to meet different demands



Tapered design for less tissue bunching and better bag closure



Exceptional safety and security



Product	Description	Packaging
PMSREB-0100	100ml, 70mm x 150mm, Single Use, Sterile	1/pk,100/bk,1000/ctn
PMSREB-0400	400ml, 120mm x 170mm, Single Use, Sterile	1/pk,100/bk,1000/ctn
PMSREB-1000	1000ml, 180mm x 220mm, Single Use, Sterile	1/pk,100/bk,1000/ctn



Disposable Wound Protector




Extra Protection, Extra Atraumatic Retraction, Ultimate Versatility


Disposable Wound Protector

Disposable wound protector is used for soft tissue and thoracic retraction, facilitates specimen removal and instruments operation. It provides 360° atraumatic retraction and reduces superficial surgical site infection following surgeries, distributes force evenly, eliminating point trauma and associated pain.

Features and Benefits

- Adapt invasive size to facilitates rapid and effortless setup.
- Provides all-around circumferential, atraumatic retraction.
- Allows visualization of wound margins.
- Provides all-around wound protection, reduces superficial surgical site infection.
- Maximizes exposure with a minimum incision size.
- Maintains moisture at the incision site.

	Model	Specification	Packaging
Regid Ring			
	PMSSG8801	Rigid Ring, Small, Sheath Length 180mm, Incision Range 25-60mm	1/pk, 10/bx, 80/ctn
	PMSSG8802	Rigid Ring, Medium, Sheath Length 180mm, Incision Range 50-90mm	
	PMSSG8803	Rigid Ring, Large, Sheath Length 250mm, Incision Range 90-140mm	
	PMSSG8804	Rigid Ring, X Large, Sheath Length 340mm, Incision Range 110-170mm	
	PMSSG8805	Rigid Ring, XX Large, Sheath Length 360mm, Incision Range 170-250mm	
Flexible Ring			
	PMSSG8831	Flexible Ring, XX-Small, Sheath Length 200mm, Incision Range 10-30mm	1/pk, 10/bx, 80/ctn
	PMSSG8832	Flexible Ring, XX-Small-Brief, Sheath Length 110mm, Incision Range 10-30mm	
	PMSSG8833	Flexible Ring, X-Small, Sheath Length 190mm, Incision Range 20-40mm	
	PMSSG8834	Flexible Ring, X-Small-Brief, Sheath Length 130mm, Incision Range 20-40mm	
	PMSSG8835	Flexible Ring, Small, Sheath Length 180mm, Incision Range 25-60mm	
	PMSSG8836	Flexible Ring, Medium, Sheath Length 180mm, Incision Range 50-90mm	
	PMSSG8837	Flexible Ring, Large, Sheath Length 250mm, Incision Range 90-140mm	
	PMSSG8838	Flexible Ring, X-Large, Sheath Length 340mm, Incision Range 110-170mm	
With Cap			
	PMSSG8861	w/cap, Small, Sheath Length 180mm, Incision Range 25-60mm	1/pk, 10/bx, 80/ctn
	PMSSG8862	w/cap, Medium, Sheath Length 180mm, Incision Range 50-90mm	



Disposable Skin Retractor












Visualization, Organization, Efficiency

Disposable Skin Retractor

The Disposable Retractor system provides a great anatomical visualization for multi-type surgeries. The variety kinds of hook placements and elastic stays maintain consistent retraction. With the Surgimed Retractor, Surgeons are free to perform other tasks with greater efficiency.

Features and Benefits

- Provide a clear, unimpeded access to the surgical site.
- Increased control, support and exposure for incisions.
- Construction eliminates the risk of detachment.
- Concave structure makes the hooks easy to access and easy to place.
- Steady Retractor maintains consistent tension throughout surgery.

	Model	Specification	Packaging
	PMSSG1601-R	Retractor, 25.3*11.3cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1602-R	Retractor, 22.0*11.3cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1603-R	Retractor, 32.5*18.3cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1604-R	Retractor, 14.1*14.1cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1605-R	Retractor, 18.6*8.9cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1606-R	Retractor, 28.3*18.3cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1607-R	Retractor, 28.3*25.0cm, Sterile Plastic, Single Use	1/pk, 20/ctn
	PMSSG1601-H-1	Sharp Hooks, 3mm, Sterile, Single Use	1/pk, 20/bx
	PMSSG1601-H-2	Sharp Hooks, 3mm, Sterile, Single Use	2/pk, 40/bx
	PMSSG1601-H-6	Sharp Hooks, 3mm, Sterile, Single Use	6/pk, 120/bx
	PMSSG1602-H-1	Sharp Hooks, 5mm, Sterile, Single Use	1/pk, 20/bx
	PMSSG1602-H-2	Sharp Hooks, 5mm, Sterile, Single Use	2/pk, 40/bx
	PMSSG1602-H-6	Sharp Hooks, 5mm, Sterile, Single Use	6/pk, 120/bx
	PMSSG1603-H-1	Blunt Hooks, 12mm, Sterile, Single Use	1/pk, 20/bx
	PMSSG1603-H-2	Blunt Hooks, 12mm, Sterile, Single Use	2/pk, 40/bx
	PMSSG1603-H-6	Blunt Hooks, 12mm, Sterile, Single Use	6/pk, 120/bx
	PMSSG1604-H-1	Two-Bifurcate Hooks, 12mm, Sterile, Single Use	1/pk, 20/bx
	PMSSG1604-H-2	Two-Bifurcate Hooks, 12mm, Sterile, Single Use	2/pk, 40/bx
	PMSSG1604-H-6	Two-Bifurcate Hooks, 12mm, Sterile, Single Use	6/pk, 120/bx



Disposable Endo Covers

Clean, Effective, Practical

Disposable Endo Covers

Surgitools Endo Covers is a latex-free, sterile, disposable protective covering for ENT endoscopes. The complete system provides a fast and effective method of reprocessing the endoscope and ensures an insertion tube covered with a clean cover for each procedure to against cross contamination.

Features and Benefits

- Optimal prevention against cross contamination.
- Time saving for sterilization process.
- Reduce probe down-time and save overhead cost.
- Minimizes exposure to germicides.

Model		Nasopharyngoscope / Endo Covers System / Sys	Packaging
	PMSES-01	Olympus Model/ENT-P4/GP	1/pk, 20/box, 200/carton
		Xion Model/EV-NG	1/pk, 20/box, 200/carton
	PMSES-02	Optim Medical/ENT 3.6 FL30/ENTity	1/pk, 20/box, 200/carton
		VNL-1190STK	1/pk, 20/box, 200/carton



Disposable Ultrasound Probe Cover

Clean, Smooth & Convenient

Disposable Ultrasound Probe Cover

The cover allows use of the transducer in scanning and needle guided procedures for multi-purpose of ultrasound diagnose, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer.

Features and Benefits

Unique material blend provides enhanced acoustic clarity and increased flexibility.

Conformal fit/shape to different types of transducer.

Rolled Product creates a clear view for transducer installation and gel application.

Prevent artifacts and provides natural nesting fit.



Model	Specification	Packaging
PMSSG2001	Sterile PE film 15.2cm tapered to 7.6*244cm, TPU film 14*30cm , Accordion. Folding, w/20g gel,Single Use	1/pk,20/ctn
PMSSG2002	Sterile PE film 15.2cm tapered to 7.6*244cm, TPU film 14*30cm , Accordion. Folding, w/0 gel,Single Use	1/pk,20/ctn
PMSSG2003	Sterile PE film 15.2cm tapered to 7.6*244cm, TPU film 14*30cm , Flat Folding, w/20g gel,Single Use	1/pk,20/ctn
PMSSG2004	Sterile TPU film 10*150cm , w/20g gel,Flat Folding,Single Use	1/pk,20/ctn
PMSSG2005	Sterile TPU film 8*12cm , w/20g gel,Flat Folding,Single Use	1/pk,20/ctn
PMSSG2006	SSterile TPU film 10*25cm , w/20g gel,Flat Folding,Single Use	1/pk,20/ctn



Disposable Trocar

Secure, Stable & Efficient

Disposable Trocar

Provides an extensive portfolio of access system to meet the variety needs of minimally invasive surgery. Surgitools trocar system with advantages of fully removable seals, first entry insufflation, and an array of fixation.

Features and Benefits

- Minimal Fascial Defect.
- Provides all-around circumferential, atraumatic retraction.
- Minimal Penetration with Rapid Insufflation.
- Repaid desufflation and specimen removal.
- Superior Abdominal Wall Retention.
- Clear Indication of Shield Position.



Model	Specification	Packaging
PMSLET1805	φ5, Stainless Steel Blade, Single Use, Sterile	1/pk,6/bx,24/ctn
PMSLET1805-T	φ5, Transparency Blade, Single Use, Sterile	1/pk,6/bx,24/ctn
PMSLET1810	φ10, Stainless Steel Blade, Single Use, Sterile	1/pk,6/bx,24/ctn
PMSLET1810-T	φ10, Transparency Blade, Single Use, Sterile	1/pk,6/bx,24/ctn
PMSLET1812	φ12, Stainless Steel Blade, Single Use, Sterile	1/pk,6/bx,24/ctn
PMSLET1812-T	φ12, Transparency Blade, Single Use, Sterile	1/pk,6/bx,24/ctn



Disposable Knot Pusher

Leggiere, Flexible, Secure



Disposable Knot Pusher


Disposable knot pusher is used during laparoscopic procedures to tighten surgical knots and to cut the suture. Knot pushers aid in simplify the tying of slipknots, simple and double knots, and can be used to form and manipulate the loop. The pusher offers an easy and safe specimen capture and removal during surgery.

Features and Benefits

- To prevent cross infection.
- Glass fiber provides extra resilience.
- Exceptional safety and security.
- Compact construction and easy to start.
- Available for multiple head types.



Product	Specification	Packaging
PMSLKP1901	Annular Type,5mm*480mm,single use,sterile	2/pk,20/bx,200/ctn
PMSLKP1902	Cavate Type,5mm*480mm,single use,sterile	2/pk,20/bx,200/ctn



Disposable Uterine Cannula

Adjustable, Effective, Easy operation



PARKER MEDICAL SUPPLIES

Disposable Uterine Cannula




Disposable Uterine Cannula provides both hydrotubation injection and uterine manipulation. The unique design allows a tight seal on the cervix and a distal extension for enhanced manipulation.

Features and Benefits

- Simple and effective
- Fully disposable and sterile packed ready for use
- Unique screw design allows better cervical seal preventing dye leakage/backflow
- Adjustable distal extension accommodates different sizes of uterus allowing uterine manipulation
- Dynamic locking mechanism for holding the tenaculum forceps



Product	Description	Packaging
PMSUCM-1810	Disposable Uterine Cannula/Manipulator, Straight Tip, Adjustable Cervical Screw, Single-use, Sterile	1/pk,10/bx,100/ctn
PMSUCM-1820	Disposable Uterine Cannula/Manipulator, Curved Tip, Adjustable Cervical Screw, Single-use, Sterile	1/pk,10/bx,100/ctn




Medical Membranes Processing

Pouches, Covers, Cases, Tubes Design and Manufacture Solution



Medical Pouches

Category	Discription	Packaging
Material	TPU, Nylon w/TPU layer, PVC etc.	
Specification	Customized, e.g. sawtooth mouth, flat mouth, w/layer, semicircle shape/V shape/Blade Shape etc.	Customized, e.g. PE Bag/Double PE Bag/Poly Bag
Process		



Probe Cover

Category	Discription	Packaging
Material	TPU,PVC etc.	
Specification	Customized, e.g. sawtooth mouth, flat mouth, w/layer, semicircle shape/V shape/Blade Shape etc.	Customized, e.g. PE Bag/Double PE Bag/Poly Bag
Process		



Tubes

Category	Discription	Packaging
Material	TPU etc.	
Specification	Customized, Horn Tube, Multi-channel Type	Customized
Process	Thermal Forming	



Liner/Air Pocket

Category	Discription	Packaging
Material	TPU, Nylon w/TPU layer, PVC etc.	
Specification	Customized	Customized
Process	Thermal Forming, High Frequency etc.	



PARKER MEDICAL SUPPLIES

CONTACT US



+1 (832) 343-3140



info@parkermedicalsupplies.com



www.parkermedicalsupplies.com

VISIT US:

2800 Post Oak Blvd, Ste 5600, Houston TX 77056, United States